

**产品名称: Mas1-like Rabbit Polyclonal Antibody**  
**产品货号: AP Rab13653**



## 产品概述 (Summary)

|                        |                                      |
|------------------------|--------------------------------------|
| 产品名称 (Production Name) | Mas1-like Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody           |
| 宿主 (Host)              | Rabbit                               |
| 应用 (Application)       | ICC/IF,ELISA                         |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                      |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | MAS1L                                                                                                  |
| 别名 (Alternative Names) | MAS1L; MRG; Mas-related G-protein coupled receptor MRG; MAS-R; MAS1-like                               |
| 基因 ID (Gene ID)        | 116511.0                                                                                               |
| 蛋白 ID (SwissProt ID)   | P35410.The antiserum was produced against synthesized peptide derived from human MAS1L. AA range:31-80 |

## 产品应用 (Application)

|                          |                                          |
|--------------------------|------------------------------------------|
| 稀释比 (Dilution Ratio)     | ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000 |
| 蛋白分子量 (Molecular Weight) |                                          |

**产品名称: Mas1-like Rabbit Polyclonal Antibody**  
**产品货号: APRab13653**




## 研究背景 (Background)

similarity:Belongs to the G-protein coupled receptor 1 family.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,similarity:Belongs to the G-protein coupled receptor 1 family.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,

## 研究领域 (Research Area)

Neuroscience

## 图片 (Image Data)



Immunofluorescence analysis of HeLa cells, using MAS1L Antibody. The picture on the right is blocked with the synthesized peptide.

## 注意事项 (Note)

For research use only .