

产品名称: MCH-2R Rabbit Polyclonal Antibody
产品货号: APRab13709



产品概述 (Summary)

产品名称 (Production Name)	MCH-2R Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCHR2
别名 (Alternative Names)	MCHR2; GPR145; SLT; Melanin-concentrating hormone receptor 2; MCH receptor 2; MCH-R2; MCHR-2; G-protein coupled receptor 145; GPRv17; MCH-2R; MCH2; MCH2R
基因 ID (Gene ID)	84539.0
蛋白 ID (SwissProt ID)	Q969V1.The antiserum was produced against synthesized peptide derived from human MCHR2. AA range:291-340

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	39kDa

产品名称: MCH-2R Rabbit Polyclonal Antibody
产品货号: APRab13709




研究背景 (Background)

function: Receptor for melanin-concentrating hormone, coupled to G proteins that activate phosphoinositide hydrolysis., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Specifically expressed in the brain, with highest levels in cerebral cortex, hippocampus and amygdala. No expression detected in the cerebellum, thalamus or hypothalamus., function: Receptor for melanin-concentrating hormone, coupled to G proteins that activate phosphoinositide hydrolysis., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Specifically expressed in the brain, with highest levels in cerebral cortex, hippocampus and amygdala. No expression detected in the cerebellum, thalamus or hypothalamus.,


研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using MCHR2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using MCHR2 antibody.

产品名称: MCH-2R Rabbit Polyclonal Antibody
产品货号: APRab13709



注意事项 (Note)

For research use only .