

产品名称: MCT13 Rabbit Polyclonal Antibody
产品货号: APRab13740



产品概述 (Summary)

产品名称 (Production Name)	MCT13 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC16A13
别名 (Alternative Names)	SLC16A13; MCT13; Monocarboxylate transporter 13; MCT 13; Solute carrier family 16 member 13
基因 ID (Gene ID)	201232.0
蛋白 ID (SwissProt ID)	Q7RTY0.The antiserum was produced against synthesized peptide derived from human MOT13. AA range:377-426

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: MCT13 Rabbit Polyclonal Antibody
产品货号: APRab13740




研究背景 (Background)


function: Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates., similarity: Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13)
family., function: Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates., similarity: Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13)
family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using MOT13 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of 3T3 cells using MCT13 Polyclonal Antibody

注意事项 (Note)

For research use only .