

产品名称: MCT14 Rabbit Polyclonal Antibody
产品货号: APRab13741



产品概述 (Summary)

产品名称 (Production Name)	MCT14 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC16A14
别名 (Alternative Names)	SLC16A14; MCT14; Monocarboxylate transporter 14; MCT 14; Solute carrier family 16 member 14
基因 ID (Gene ID)	151473.0
蛋白 ID (SwissProt ID)	Q7RTX9.The antiserum was produced against synthesized peptide derived from human MOT14. AA range:131-180

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	48kDa

产品名称: MCT14 Rabbit Polyclonal Antibody
产品货号: APRab13741





研究背景 (Background)

function: Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates., similarity: Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13)
family., function: Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates., similarity: Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13)
family.,

研究领域 (Research Area)

Signal Transduction; Metabolism; Plasma Membrane; Channels

图片 (Image Data)



注意事项 (Note)

For research use only .