

产品名称: Melan-A Rabbit Polyclonal Antibody
产品货号: APRab13814



产品概述 (Summary)

产品名称 (Production Name)	Melan-A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLANA
别名 (Alternative Names)	MLANA; MART1; Melanoma antigen recognized by T-cells 1; MART-1; Antigen LB39-AA; Antigen SK29-AA; Protein Melan-A
基因 ID (Gene ID)	2315.0
蛋白 ID (SwissProt ID)	Q16655.The antiserum was produced against synthesized peptide derived from human MART-1. AA range:41-90

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: Melan-A Rabbit Polyclonal Antibody
产品货号: AP Rab13814




研究背景 (Background)

tissue specificity: Expression is restricted to melanoma and melanocyte cell lines and retina., tissue specificity: Expression is restricted to melanoma and melanocyte cell lines and retina.,


研究领域 (Research Area)

Tags & Cell Markers; Cell Type Markers; Tumor Associated; Cancer; Tumor immunology; Tumor-associated antigens

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using MART-1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from NIH/3T3 cells, using MART-1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Melan-A Rabbit Polyclonal Antibody
产品货号: APRab13814



Western Blot analysis of 3T3 cells using Melan-A Polyclonal Antibody

注意事项 (Note)

For research use only .