

产品名称: MIA Rabbit Polyclonal Antibody
产品货号: APRab13881



产品概述 (Summary)

产品名称 (Production Name)	MIA Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MIA
别名 (Alternative Names)	MIA; Melanoma-derived growth regulatory protein; Melanoma inhibitory activity protein
基因 ID (Gene ID)	8190.0
蛋白 ID (SwissProt ID)	Q16674. The antiserum was produced against synthesized peptide derived from human MIA. AA range: 82-131

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: MIA Rabbit Polyclonal Antibody
产品货号: APRab13881


研究背景 (Background)

function: Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas., PTM: May possess two intramolecular disulfide bonds., similarity: Belongs to the MIA/OTOR family., similarity: Contains 1 SH3 domain., tissue specificity: All malignant melanoma cell lines tested and infrequently in glioma cell lines., function: Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas., PTM: May possess two intramolecular disulfide bonds., similarity: Belongs to the MIA/OTOR family., similarity: Contains 1 SH3 domain., tissue specificity: All malignant melanoma cell lines tested and infrequently in glioma cell lines.,

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human heart tissue, using MIA Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .