

产品名称: MIG Rabbit Polyclonal Antibody
产品货号: AP Rab13902



产品概述 (Summary)

产品名称 (Production Name)	MIG Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CXCL9 CMK MIG SCYB9
别名 (Alternative Names)	C-X-C motif chemokine 9 (Gamma-interferon-induced monokine; Monokine induced by interferon-gamma; HuMIG; MIG; Small-inducible cytokine B9)
基因 ID (Gene ID)	4283.0
蛋白 ID (SwissProt ID)	Q07325.Synthetic peptide from human protein at AA range: 90-125

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: MIG Rabbit Polyclonal Antibody
产品货号: APRab13902

This antimicrobial gene encodes a protein thought to be involved in T cell trafficking. The encoded protein binds to C-X-C motif chemokine 3 and is a chemoattractant for lymphocytes but not for neutrophils. [provided by RefSeq, Sep 2014],function:Cytokine that affects the growth, movement, or activation state of cells that participate in immune and inflammatory response. Chemotactic for activated T-cells. Binds to CXCR3.,induction:By interferon gamma. The induction is enhanced by TNF-alpha in dermal fibroblasts and vein endothelial cells.,online information:CXCL9 entry,similarity:Belongs to the intercrine alpha (chemokine CxC) family.,


研究领域 (Research Area)

Cytokine-cytokine receptor interaction;Chemokine;Toll_Like;

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

注意事项 (Note)

For research use only .