

产品名称: MIP-T3 Rabbit Polyclonal Antibody
产品货号: AP Rab13914



产品概述 (Summary)

产品名称 (Production Name)	MIP-T3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAF3IP1 TRAF3IP1; MIPT3; TRAF3-interacting protein 1; Interleukin-13 receptor alpha
别名 (Alternative Names)	1-binding protein 1; Microtubule-interacting protein associated with TRAF3; MIP-T3
基因 ID (Gene ID)	26146.0
蛋白 ID (SwissProt ID)	Q8TDR0.The antiserum was produced against synthesized peptide derived from human MIPT3. AA range:221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	78kDa

产品名称: MIP-T3 Rabbit Polyclonal Antibody
产品货号: APRab13914




研究背景 (Background)

function: Play an inhibitory role on IL13 signaling by binding to IL13RA1. Involved in suppression of IL13-induced STAT6 phosphorylation, transcriptional activity and DNA-binding. Recruits TRAF3 and DISC1 to the microtubules., similarity: Belongs to the TRAF3IP1 family., subcellular location: Microtubules., subunit: Interacts with IL13RA1. Binds to microtubules, TRAF3 and DISC1., tissue specificity: Ubiquitous., function: Play an inhibitory role on IL13 signaling by binding to IL13RA1. Involved in suppression of IL13-induced STAT6 phosphorylation, transcriptional activity and DNA-binding. Recruits TRAF3 and DISC1 to the microtubules., similarity: Belongs to the TRAF3IP1 family., subcellular location: Microtubules., subunit: Interacts with IL13RA1. Binds to microtubules, TRAF3 and DISC1., tissue specificity: Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of MCF7 cells, using MIP-T3 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .