

产品名称: MN1 Rabbit Polyclonal Antibody
产品货号: APRab14006



产品概述 (Summary)

产品名称 (Production Name)	MN1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MN1
别名 (Alternative Names)	MN1; Probable tumor suppressor protein MN1
基因 ID (Gene ID)	4330.0
蛋白 ID (SwissProt ID)	Q10571.The antiserum was produced against synthesized peptide derived from human MN1. AA range:821-870

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	135kDa

研究背景 (Background)


产品名称: MN1 Rabbit Polyclonal Antibody
产品货号: APRab14006




Meningioma 1 (MN1) contains two sets of CAG repeats. It is disrupted by a balanced translocation (4;22) in a meningioma, and its inactivation may contribute to meningioma pathogenesis. [provided by RefSeq, Jul 2008],disease:A chromosomal aberration involving MN1 may be a cause of acute myeloid leukemia (AML). Translocation t(12;22)(p13;q11) with TEL,disease:Defects in MN1 may be a cause of meningiomas, slowly growing benign tumors derived from the arachnoidal cap cells of the leptomeninges, the soft coverings of the brain and spinal cord. Meningiomas are believed to be the most common primary tumors of the central nervous system in man.,function:May play a role in tumor suppression.,tissue specificity:Ubiquitously expressed. Highest levels in skeletal muscle.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COLO cells, using MN1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MN1 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .