

产品名称: Mob1A Rabbit Polyclonal Antibody
产品货号: APRab14013



产品概述 (Summary)

产品名称 (Production Name)	Mob1A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MOB1B
别名 (Alternative Names)	MOB1B; MOB4A; MOBKL1A; MOB kinase activator 1B; Mob1 homolog 1A; Mob1A; Mob1B; Mps one binder kinase activator-like 1A
基因 ID (Gene ID)	92597.0
蛋白 ID (SwissProt ID)	Q7L9L4.The antiserum was produced against synthesized peptide derived from human MOL1A. AA range:101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	25kDa

产品名称: Mob1A Rabbit Polyclonal Antibody
产品货号: APRab14013




研究背景 (Background)


The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],function:Stimulates the kinase activity of STK38L, similarity:Belongs to the MOB1/phocein family, subunit: Binds STK38L,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from 293, COS7, and HUVEC cells, treated with IFN 2500U/ml 30', using MOL1A Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Mob1A Polyclonal Antibody

注意事项 (Note)

For research use only .