

产品名称: Mob3A/B Rabbit Polyclonal Antibody
产品货号: AP Rab14015



产品概述 (Summary)

产品名称 (Production Name)	Mob3A/B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MMOBKL2A/MMOBKL2B MOB3A; MOBKL2A; MOB kinase activator 3A; MOB-LAK; Mob1 homolog 2A;
别名 (Alternative Names)	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B
基因 ID (Gene ID)	126308.0
蛋白 ID (SwissProt ID)	Q96BX8.Synthesized peptide derived from Mob3A/B . at AA range: 1-80

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	25kDa

产品名称: Mob3A/B Rabbit Polyclonal Antibody
产品货号: AP Rab14015




研究背景 (Background)


function: May regulate the activity of kinases., sequence caution: Chimera. There is an inclusion of DNA from another gene in position 141., similarity: Belongs to the MOB1/phocein family., function: May regulate the activity of kinases., sequence caution: Chimera. There is an inclusion of DNA from another gene in position 141., similarity: Belongs to the MOB1/phocein family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 3T3 using Mob3A/B antibody.. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200

注意事项 (Note)

For research use only .