

产品名称: MOV10L1 Rabbit Polyclonal Antibody
产品货号: APRab14048



产品概述 (Summary)

产品名称 (Production Name)	MOV10L1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MOV10L1
别名 (Alternative Names)	MOV10L1; Putative helicase Mov10l1; Moloney leukemia virus 10-like protein 1; MOV10-like protein 1
基因 ID (Gene ID)	54456.0
蛋白 ID (SwissProt ID)	Q9BXT6.The antiserum was produced against synthesized peptide derived from human MOV10L1. AA range:318-367

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	140kDa

产品名称: MOV10L1 Rabbit Polyclonal Antibody
产品货号: AP Rab14048




研究背景 (Background)


Mov10 RISC complex RNA helicase like 1(MOV10L1) Homo sapiens This gene is similar to a mouse gene that encodes a putative RNA helicase and shows testis-specific expression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009], alternative products: Additional isoforms seem to exist, function: Putative RNA helicase. Isoform 1 may play a role in male germ cell development., similarity: Belongs to the DNA2/NAM7 helicase family. SDE3 subfamily., tissue specificity: Isoform 1 is specifically expressed in testis.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa and RAW264.7 cells, using MOV10L1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using MOV10L1 antibody.

产品名称: MOV10L1 Rabbit Polyclonal Antibody
产品货号: APRab14048



Western Blot analysis of various cells using MOV10L1 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .