

产品名称: MPP9 Rabbit Polyclonal Antibody
产品货号: APRab14064



产品概述 (Summary)

产品名称 (Production Name)	MPP9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MPHOSPH9
别名 (Alternative Names)	MPHOSPH9; MPP9; M-phase phosphoprotein 9
基因 ID (Gene ID)	10198.0
蛋白 ID (SwissProt ID)	Q99550.The antiserum was produced against synthesized peptide derived from human MPHOSPH9. AA range:539-588

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	78kDa

研究背景 (Background)


产品名称: MPP9 Rabbit Polyclonal Antibody
产品货号: APRab14064




PTM:Phosphorylated in M (mitotic) phase.,PTM:Phosphorylated in M (mitotic) phase.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of MPHOSPH9 Antibody. The lane on the right is blocked with the MPHOSPH9 peptide.



Western blot analysis of the lysates from COLO205 cells using MPHOSPH9 antibody.



Western Blot analysis of various cells using MPP9 Polyclonal Antibody

产品名称: MPP9 Rabbit Polyclonal Antibody
产品货号: APRab14064



注意事项 (Note)

For research use only .