

产品名称: mPR α Rabbit Polyclonal Antibody
产品货号: AP Rab14069



产品概述 (Summary)

产品名称 (Production Name)	mPR α Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAQR7 PAQR7; MRPA; Membrane progesterin receptor alpha; mPR alpha; Progesterin and adipoQ receptor family member 7; Progesterin and adipoQ receptor family member VII
别名 (Alternative Names)	
基因 ID (Gene ID)	164091.0
蛋白 ID (SwissProt ID)	Q86WK9.The antiserum was produced against synthesized peptide derived from human MPRA. AA range:297-346

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	40kDa

产品名称: mPR α Rabbit Polyclonal Antibody
产品货号: AP Rab14069




研究背景 (Background)


function: Steroid membrane receptor. Binds progesterone in vitro. May be involved in oocyte maturation., similarity: Belongs to the ADIPOR family., tissue specificity: Expressed in ovary, testis, placenta, uterus and bladder., function: Steroid membrane receptor. Binds progesterone in vitro. May be involved in oocyte maturation., similarity: Belongs to the ADIPOR family., tissue specificity: Expressed in ovary, testis, placenta, uterus and bladder.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from 293 and COLO cells, using MPRA Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of 293 cells using mPR α Polyclonal Antibody

注意事项 (Note)

For research use only .