

产品名称: mPR γ Rabbit Polyclonal Antibody
产品货号: AP Rab14070



产品概述 (Summary)

产品名称 (Production Name)	mPR γ Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAQR5 PAQR5; MPRG; Membrane progesterin receptor gamma; mPR gamma; Progesterin and adipoQ receptor family member 5; Progesterin and adipoQ receptor family member V
别名 (Alternative Names)	
基因 ID (Gene ID)	54852.0
蛋白 ID (SwissProt ID)	Q9NXX6.The antiserum was produced against synthesized peptide derived from human MPRG. AA range:281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: mPRy Rabbit Polyclonal Antibody
产品货号: APRab14070




研究背景 (Background)

function: Steroid membrane receptor. Binds progesterone. May be involved in oocyte maturation., similarity: Belongs to the ADIPOR family., tissue specificity: Expressed in the kidney, colon, adrenal and lung., function: Steroid membrane receptor. Binds progesterone. May be involved in oocyte maturation., similarity: Belongs to the ADIPOR family., tissue specificity: Expressed in the kidney, colon, adrenal and lung.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using MPRG Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .