

产品名称: MRGF Rabbit Polyclonal Antibody
产品货号: APRab14087



产品概述 (Summary)

产品名称 (Production Name)	MRGF Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRGPRF MRGPRF; GPR140; GPR168; MRGF; PSEC0142; Mas-related G-protein coupled receptor member F; Mas-related gene F protein; G-protein coupled receptor 140; G-protein coupled receptor 168
别名 (Alternative Names)	
基因 ID (Gene ID)	116535.0
蛋白 ID (SwissProt ID)	Q96AM1.The antiserum was produced against synthesized peptide derived from human MRGRF. AA range:271-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	40kDa

产品名称: MRGF Rabbit Polyclonal Antibody
产品货号: APRab14087




研究背景 (Background)

function:Orphan receptor. May bind to a neuropeptide and may regulate nociceptor function and/or development, including the sensation or modulation of pain.,similarity:Belongs to the G-protein coupled receptor 1


family.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,function:Orphan receptor. May bind to a neuropeptide and may regulate nociceptor function and/or development, including the sensation or modulation of pain.,similarity:Belongs to the G-protein coupled receptor 1 family.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using MRGRF Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MRGF Polyclonal Antibody

注意事项 (Note)

For research use only .

产品名称: MRGF Rabbit Polyclonal Antibody
产品货号: APRab14087

