

产品名称: MRGG Rabbit Polyclonal Antibody
产品货号: APRab14088



产品概述 (Summary)

产品名称 (Production Name)	MRGG Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRGPRG
别名 (Alternative Names)	MRGPRG; GPR169; MRGG; Mas-related G-protein coupled receptor member G; G-protein coupled receptor 169
基因 ID (Gene ID)	386746.0
蛋白 ID (SwissProt ID)	Q86SM5.The antiserum was produced against synthesized peptide derived from human MRGRG. AA range:231-280

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

产品名称: MRGG Rabbit Polyclonal Antibody
产品货号: APRab14088




研究背景 (Background)


function:Orphan receptor. May regulate nociceptor function and/or development, including the sensation or modulation of pain.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,function:Orphan receptor. May regulate nociceptor function and/or development, including the sensation or modulation of pain.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using MRGRG Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from HUVEC cells, using MRGRG Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: MRGG Rabbit Polyclonal Antibody
产品货号: APRab14088



Western blot analysis of the lysates from HepG2 cells using MRGRG antibody.



Western Blot analysis of various cells using MRGG Polyclonal Antibody

注意事项 (Note)

For research use only .