

产品名称: MRP-L11 Rabbit Polyclonal Antibody
产品货号: APRab14109



产品概述 (Summary)

产品名称 (Production Name)	MRP-L11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL11
别名 (Alternative Names)	MRPL11; CGI-113; 39S ribosomal protein L11; mitochondrial; L11mt; MRP-L11
基因 ID (Gene ID)	65003.0
蛋白 ID (SwissProt ID)	Q9Y3B7.The antiserum was produced against synthesized peptide derived from human MRPL11. AA range:21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	21kDa

研究背景 (Background)

产品名称: MRP-L11 Rabbit Polyclonal Antibody
产品货号: AP Rab14109




This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12. [provided by RefSeq, May 2014],similarity:Belongs to the ribosomal protein L11P family.,


研究领域 (Research Area)

Cell Biology

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using MRPL11 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human cervix carcinoma tissue, using MRPL11 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: MRP-L11 Rabbit Polyclonal Antibody
产品货号: APRab14109



Western blot analysis of lysates from COS cells, using MRPL11 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MRP-L11 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .