

产品名称: MRP-L14 Rabbit Polyclonal Antibody
产品货号: APRab14112



产品概述 (Summary)

产品名称 (Production Name)	MRP-L14 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL14
别名 (Alternative Names)	MRPL14; MRPL32; RPML32; 39S ribosomal protein L14; mitochondrial; L14mt; MRP-L14; 39S ribosomal protein L32, mitochondrial; L32mt; MRP-L32
基因 ID (Gene ID)	64928.0
蛋白 ID (SwissProt ID)	Q6P1L8.The antiserum was produced against synthesized peptide derived from human MRPL14. AA range:96-145

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	16kDa

产品名称: MRP-L14 Rabbit Polyclonal Antibody
产品货号: APRab14112




研究背景 (Background)


This nuclear gene encodes a protein component of the 39S subunit of the mitochondrial ribosome. A pseudogene of this gene is found on chromosome 17. [provided by RefSeq, Jan 2016],similarity:Belongs to the ribosomal protein L14P family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from 293, HeLa, and COLO cells, using MRPL14 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using MRPL14 antibody.



产品名称: MRP-L14 Rabbit Polyclonal Antibody
产品货号: APRab14112



Western Blot analysis of various cells using MRP-L14 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .