

产品名称: MRP-L18 Rabbit Polyclonal Antibody
产品货号: APRab14116



产品概述 (Summary)

产品名称 (Production Name)	MRP-L18 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL18
别名 (Alternative Names)	MRPL18; HSPC071; 39S ribosomal protein L18; mitochondrial; L18mt; MRP-L18
基因 ID (Gene ID)	29074.0
蛋白 ID (SwissProt ID)	Q9H0U6.The antiserum was produced against synthesized peptide derived from human MRPL18. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	21kDa

产品名称: MRP-L18 Rabbit Polyclonal Antibody
产品货号: APRab14116




研究背景 (Background)


This nuclear gene encodes a protein component of the larger 39S subunit of mitochondrial ribosome. This protein may also aid in the import of nuclear-encoded 5S rRNA into mitochondria. Alternative splicing results in multiple transcript variants, most of which are not predicted to encode a protein. A pseudogene of this gene is found on chromosome 16. [provided by RefSeq, Jan 2016],similarity:Belongs to the ribosomal protein L18P family.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HT-29 cells, using MRPL18 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using MRPL18 antibody.

产品名称: MRP-L18 Rabbit Polyclonal Antibody
产品货号: APRab14116



Western Blot analysis of various cells using MRP-L18 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .