

产品名称: MRP-L44 Rabbit Polyclonal Antibody
产品货号: APRab14130



产品概述 (Summary)

产品名称 (Production Name)	MRP-L44 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL44
别名 (Alternative Names)	MRPL44; 39S ribosomal protein L44; mitochondrial; L44mt; MRP-L44
基因 ID (Gene ID)	65080.0
蛋白 ID (SwissProt ID)	Q9H9J2.The antiserum was produced against synthesized peptide derived from human MRPL44. AA range:221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)


产品名称: MRP-L44 Rabbit Polyclonal Antibody
产品货号: APRab14130




Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. [provided by RefSeq, Jul 2008],function:Component of the 39S subunit of mitochondrial ribosome.,similarity:Contains 1 DRBM (double-stranded RNA-binding) domain.,similarity:Contains 1 RNase III domain.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of HUVEC cells, using MRPL44 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human heart tissue, using MRPL44 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: MRP-L44 Rabbit Polyclonal Antibody
产品货号: APRab14130



Western blot analysis of lysates from COS7 cells, using MRPL44 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MRP-L44 Polyclonal Antibody

注意事项 (Note)

For research use only .