

产品名称: MRP-L9 Rabbit Polyclonal Antibody
产品货号: APRab14139



产品概述 (Summary)

产品名称 (Production Name)	MRP-L9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL9
别名 (Alternative Names)	MRPL9; 39S ribosomal protein L9; mitochondrial; L9mt; MRP-L9
基因 ID (Gene ID)	65005.0
蛋白 ID (SwissProt ID)	Q9BYD2.The antiserum was produced against synthesized peptide derived from human MRPL9. AA range:211-260

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	30kDa

研究背景 (Background)


产品名称: MRP-L9 Rabbit Polyclonal Antibody
产品货号: AP Rab14139




This is a nuclear gene encoding a protein component of the 39S subunit of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 8. [provided by RefSeq, Jul 2014],similarity:Belongs to the ribosomal protein L9P family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COLO cells, using MRPL9 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using MRPL9 antibody.



产品名称: MRP-L9 Rabbit Polyclonal Antibody
产品货号: APRab14139



Western Blot analysis of various cells using MRP-L9 Polyclonal Antibody

注意事项 (Note)

For research use only .