

产品名称: MSS4 Rabbit Polyclonal Antibody
产品货号: APRab14185



产品概述 (Summary)

产品名称 (Production Name)	MSS4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RABIF MSS4 RASGRF3
别名 (Alternative Names)	
基因 ID (Gene ID)	5877.0
蛋白 ID (SwissProt ID)	P47224.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	13kDa

研究背景 (Background)

This gene encodes a member of the SCE4/YPT1/RAB family of small GTP-binding proteins that are involved in the regulation of intracellular vesicular transport. This protein stimulates GTP-GDP exchange in SEC4, and to a lesser extent in


产品名称: MSS4 Rabbit Polyclonal Antibody
产品货号: APRab14185



YPT1 and RAB3A, and may play a general role in vesicular transport. [provided by RefSeq, Oct 2011],function:Guanine-nucleotide-releasing protein that acts on members of the SCE4/YPT1/RAB subfamily. Stimulates GDP release from both YPT1 and RAB3A, but is less active on these proteins than on the SEC4 protein. Might play a general role in vesicular transport.,similarity:Belongs to the DSS4/MSS4 family.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .