

产品名称: MUM1 Rabbit Polyclonal Antibody
产品货号: APRab14243



产品概述 (Summary)

产品名称 (Production Name)	MUM1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MUM1 EXPAND1
别名 (Alternative Names)	melanoma associated antigen (mutated) 1
基因 ID (Gene ID)	84939.0
蛋白 ID (SwissProt ID)	Q2TAK8.The antiserum was produced against synthesized peptide derived from the C-terminal region of human MUM1. AA range:661-710

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	78kDa

研究背景 (Background)


产品名称: MUM1 Rabbit Polyclonal Antibody
产品货号: APRab14243



PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the MUM1 family.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the MUM1 family.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of mouse-kidney mouse-heart mouse-brain cells using MUM1 Polyclonal Antibody diluted at 1:1000.
Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .