

产品名称: MxB Rabbit Polyclonal Antibody
产品货号: AP Rab14254



产品概述 (Summary)

产品名称 (Production Name)	MxB Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MX2 MX2; Interferon-induced GTP-binding protein Mx2; Interferon-regulated resistance GTP-binding protein MxB; Myxovirus resistance protein 2; p78-related protein
别名 (Alternative Names)	
基因 ID (Gene ID)	4600.0
蛋白 ID (SwissProt ID)	P20592.Synthesized peptide derived from the Internal region of human MxB.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	82kDa

产品名称: MxB Rabbit Polyclonal Antibody
产品货号: APRab14254




研究背景 (Background)

The protein encoded by this gene has a nuclear and a cytoplasmic form and is a member of both the dynamin family and the family of large GTPases. The nuclear form is localized in a granular pattern in the heterochromatin region beneath the nuclear envelope. A nuclear localization signal (NLS) is present at the amino terminal end of the nuclear form but is lacking in the cytoplasmic form due to use of an alternate translation start codon. This protein is upregulated by interferon-alpha but does not contain the antiviral activity of a similar myxovirus resistance protein 1. [provided by RefSeq, Jul 2008],induction:By interferons.,similarity:Belongs to the dynamin family.,similarity:Contains 1 GED domain.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of 293 cells using MxB Polyclonal Antibody

注意事项 (Note)

For research use only .