

产品名称: MXRA5 Rabbit Polyclonal Antibody
产品货号: APRab14256



产品概述 (Summary)

产品名称 (Production Name)	MXRA5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MXRA5
别名 (Alternative Names)	MXRA5
基因 ID (Gene ID)	25878.0
蛋白 ID (SwissProt ID)	Q9NR99.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:300, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	311kDa

研究背景 (Background)


This gene encodes one of the matrix-remodelling associated proteins. This protein contains 7 leucine-rich repeats and 12 immunoglobulin-like C2-type domains related to perlecan. This gene has a pseudogene on chromosome Y. [provided by

产品名称: MXRA5 Rabbit Polyclonal Antibody
产品货号: APRab14256

RefSeq, Mar 2010],developmental stage:Over-expressed in centenarians. Expression is reduced from young to old but increased from old to centenarians.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 12 Ig-like C2-type (immunoglobulin-like) domains.,similarity:Contains 7 LRR (leucine-rich) repeats.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .