

产品名称: MYCNOS Rabbit Polyclonal Antibody
产品货号: APRab14269



产品概述 (Summary)

产品名称 (Production Name)	MYCNOS Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYCNOS CYMN NCYM
别名 (Alternative Names)	N-cym protein (N-myc opposite strand)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P40205.Synthesized peptide derived from human MYCNOS AA range: 2-52

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

This gene is transcribed in antisense to the v-myc avian myelocytomatosis viral oncogene neuroblastoma derived homolog


产品名称: MYCNOS Rabbit Polyclonal Antibody
产品货号: APRab14269



gene (MYCN). It is thought to encode a small, novel protein that stabilizes MYCN, prevents apoptosis, and promotes cell proliferation. Transcripts at this locus may also act directly as functional RNAs to recruit transcriptional regulators to the promoter of MYCN and stimulate transcription of this oncogene. This gene therefore functions through both RNA and protein products. [provided by RefSeq, Aug 2016],developmental stage:Expressed during fetal development, as well as in tumor cell lines containing amplified N-myc loci, where it is expressed at very high levels.,function:May have a functional role during normal fetal development.,tissue specificity:Expressed in the fetal brain, lung, liver and kidney at varying low levels.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200 (4° overnight) .
2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

注意事项 (Note)

For research use only .