

产品名称: NACAD Rabbit Polyclonal Antibody
产品货号: APRab14390



产品概述 (Summary)

产品名称 (Production Name)	NACAD Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NACAD
别名 (Alternative Names)	NACAD; KIAA0363; NAC-alpha domain-containing protein 1
基因 ID (Gene ID)	23148.0
蛋白 ID (SwissProt ID)	O15069.The antiserum was produced against synthesized peptide derived from human NACAD. AA range:376-425

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	140kDa

研究背景 (Background)


产品名称: NACAD Rabbit Polyclonal Antibody
产品货号: APRab14390




function: May prevent inappropriate targeting of non-secretory polypeptides to the endoplasmic reticulum (ER). May bind to nascent polypeptide chains as they emerge from the ribosome and block their interaction with the signal recognition particle (SRP), which normally targets nascent secretory peptides to the ER. May also reduce the inherent affinity of ribosomes for protein translocation sites in the ER membrane (M sites)., similarity: Belongs to the NAC-alpha family., similarity: Contains 1 NAC-A/B (NAC-alpha/beta) domain., function: May prevent inappropriate targeting of non-secretory polypeptides to the endoplasmic reticulum (ER). May bind to nascent polypeptide chains as they emerge from the ribosome and block their interaction with the signal recognition particle (SRP), which normally targets nascent secretory peptides to the ER. May also reduce the inherent affinity of ribosomes for protein translocation sites in the ER membrane (M sites)., similarity: Belongs to the NAC-alpha family., similarity: Contains 1 NAC-A/B (NAC-alpha/beta) domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from Jurkat cells, using NACAD Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using NACAD antibody.

产品名称: NACAD Rabbit Polyclonal Antibody
产品货号: APRab14390



Western Blot analysis of Jurkat cells using NACAD Polyclonal Antibody

注意事项 (Note)

For research use only .