

产品名称: NBAS Rabbit Polyclonal Antibody
产品货号: APRab14421



产品概述 (Summary)

产品名称 (Production Name)	NBAS Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NBAS
别名 (Alternative Names)	NAG
基因 ID (Gene ID)	51594.0
蛋白 ID (SwissProt ID)	A2RRP1.Synthesized peptide derived from human protein . at AA range: 1880-1960

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:300, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	260kDa

研究背景 (Background)

This gene encodes a protein with two leucine zipper domains, a ribosomal protein S14 signature domain and a Sec39 like

