

产品名称: NCAM2 Rabbit Polyclonal Antibody
产品货号: APRab14437



产品概述 (Summary)

产品名称 (Production Name)	NCAM2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCAM2
别名 (Alternative Names)	NCAM2; NCAM21; Neural cell adhesion molecule 2; N-CAM-2; NCAM-2
基因 ID (Gene ID)	4685.0
蛋白 ID (SwissProt ID)	O15394.The antiserum was produced against synthesized peptide derived from human NCAM2. AA range:264-313

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	100kDa

研究背景 (Background)

产品名称: NCAM2 Rabbit Polyclonal Antibody
产品货号: APRab14437




The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein and may function in selective fasciculation and zone-to-zone projection of the primary olfactory axons. [provided by RefSeq, Jul 2008],function:May play important roles in selective fasciculation and zone-to-zone projection of the primary olfactory axons.,similarity:Contains 2 fibronectin type-III domains.,similarity:Contains 5 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Expressed most strongly in adult and fetal brain.,

研究领域 (Research Area)

Cell adhesion molecules (CAMs);Prion diseases;

图片 (Image Data)



Western blot analysis of lysates from MCF-7 and HeLa cells, using NCAM2 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .