

产品名称: NDF4 Rabbit Polyclonal Antibody
产品货号: AP Rab14467



产品概述 (Summary)

产品名称 (Production Name)	NDF4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEUROD4 ATH3 ATOH3 BHLHA4
别名 (Alternative Names)	
基因 ID (Gene ID)	58158.0
蛋白 ID (SwissProt ID)	Q9HD90.Synthesized peptide derived from human protein . at AA range: 30-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

function:Appears to mediate neuronal differentiation.,PTM:Serine or threonine phosphorylation within the basic region may


产品名称: NDF4 Rabbit Polyclonal Antibody
产品货号: APRab14467



regulate neurogenic activity.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,function:Appears to mediate neuronal differentiation.,PTM:Serine or threonine phosphorylation within the basic region may regulate neurogenic activity.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .