

产品名称: NDF6 Rabbit Polyclonal Antibody
产品货号: APRab14468



产品概述 (Summary)

产品名称 (Production Name)	NDF6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEUROD6 ATOH2 BHLHA2 My051
别名 (Alternative Names)	
基因 ID (Gene ID)	63974.0
蛋白 ID (SwissProt ID)	Q96NK8.Synthesized peptide derived from human protein . at AA range: 230-310

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)

This gene is a member of the NEUROD family of basic helix-loop-helix transcription factors. The encoded protein may be

产品名称: NDF6 Rabbit Polyclonal Antibody
产品货号: APRab14468




involved in the development and differentiation of the nervous system. [provided by RefSeq, Nov 2012],function:Activates E box-dependent transcription in collaboration with TCF3/E47. May be a trans-acting factor involved in the development and maintenance of the mammalian nervous system. Transactivates the promoter of its own gene.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Domain Families; HLH / Leucine Zipper; HLH; Neuroscience; Neurology process; Neurogenesis; Transcription Factors; Signal Transduction; Antibodies; New products Recombinant

图片 (Image Data)



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

注意事项 (Note)

For research use only .