

产品名称: NDRG3 Rabbit Polyclonal Antibody
产品货号: AP Rab14479



产品概述 (Summary)

产品名称 (Production Name)	NDRG3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDRG3
别名 (Alternative Names)	NDRG3; Protein NDRG3; N-myc downstream-regulated gene 3 protein
基因 ID (Gene ID)	57446.0
蛋白 ID (SwissProt ID)	Q9UGV2.The antiserum was produced against synthesized peptide derived from human NDRG3. AA range:206-255

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

产品名称: NDRG3 Rabbit Polyclonal Antibody
产品货号: AP Rab14479




similarity:Belongs to the NDRG family.,tissue specificity:Ubiquitous. Highly expressed in brain.,similarity:Belongs to the NDRG family.,tissue specificity:Ubiquitous. Highly expressed in brain.,


研究领域 (Research Area)

Cell Biology; Cell Cycle; Cell differentiation; Cancer; Tumor biomarkers

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using NDRG3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 cells, using NDRG3 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .