

产品名称: NDUFA4L2 Rabbit Polyclonal Antibody
产品货号: APRab14499



产品概述 (Summary)

产品名称 (Production Name)	NDUFA4L2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFA4L2
别名 (Alternative Names)	NDUFA4L2; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2; NADH-ubiquinone oxidoreductase MLRQ subunit homolog; NUOMS
基因 ID (Gene ID)	56901.0
蛋白 ID (SwissProt ID)	Q9NRX3.The antiserum was produced against synthesized peptide derived from human NDUFA4L2. AA range:38-87

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: NDUFA4L2 Rabbit Polyclonal Antibody
产品货号: APRab14499




研究背景 (Background)

similarity:Belongs to the complex I NDUFA4 subunit family.,similarity:Belongs to the complex I NDUFA4 subunit family.,


研究领域 (Research Area)

Oxidative phosphorylation;Alzheimer's disease;Parkinson's disease;Huntington's disease;

图片 (Image Data)



Western blot analysis of lysates from HepG2, HUVEC, COLO, and Jurkat cells, using NDUFA4L2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using NDUFA4L2 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .