

产品名称: NDUFS3 Rabbit Polyclonal Antibody
产品货号: APRab14516



产品概述 (Summary)

产品名称 (Production Name)	NDUFS3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFS3
别名 (Alternative Names)	NDUFS3; NADH dehydrogenase [ubiquinone] iron-sulfur protein 3; mitochondrial; Complex I-30kD; CI-30kD; NADH-ubiquinone oxidoreductase 30 kDa subunit
基因 ID (Gene ID)	4722.0
蛋白 ID (SwissProt ID)	O75489.The antiserum was produced against synthesized peptide derived from human NDUFS3. AA range:117-166

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: NDUFS3 Rabbit Polyclonal Antibody
产品货号: APRab14516




研究背景 (Background)

This gene encodes one of the iron-sulfur protein (IP) components of mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mutations in this gene are associated with Leigh syndrome resulting from mitochondrial complex I deficiency. [provided by RefSeq, Apr 2009],catalytic activity:NADH + acceptor = NAD(+) + reduced acceptor.,catalytic activity:NADH + ubiquinone = NAD(+) + ubiquinol.,function:Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.,similarity:Belongs to the complex I 30 kDa subunit family.,subunit:Mammalian complex I is composed of 45 different subunits.,


研究领域 (Research Area)

Oxidative phosphorylation;Alzheimer's disease;Parkinson's disease;Huntington's disease;

图片 (Image Data)




Western blot analysis of lysates from mouse brain, using NDUFS3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using NDUFS3 antibody.

产品名称: NDUFS3 Rabbit Polyclonal Antibody
产品货号: APRab14516



Western Blot analysis of 293T cells using NDUFS3 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .