

产品名称: NEEP21 Rabbit Polyclonal Antibody
产品货号: APRab14541



产品概述 (Summary)

产品名称 (Production Name)	NEEP21 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NSG1
别名 (Alternative Names)	NSG1; D4S234; Neuron-specific protein family member 1; Brain neuron cytoplasmic protein 1
基因 ID (Gene ID)	27065.0
蛋白 ID (SwissProt ID)	P42857.The antiserum was produced against synthesized peptide derived from human NSG1. AA range:112-161

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	26kDa

产品名称: NEEP21 Rabbit Polyclonal Antibody
产品货号: APRab14541




研究背景 (Background)


similarity:Belongs to the NSG family.,similarity:Belongs to the NSG family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COLO cells, using NSG1 Antibody. The lane on the right is blocked with the synthesized peptide.




Western blot analysis of the lysates from K562 cells using NSG1 antibody.



Western Blot analysis of various cells using NEEP21 Polyclonal Antibody diluted at 1: 2000

产品名称: NEEP21 Rabbit Polyclonal Antibody
产品货号: APRab14541



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .