

产品名称: NEGR1 Rabbit Polyclonal Antibody
产品货号: APRab14544



产品概述 (Summary)

产品名称 (Production Name)	NEGR1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEGR1
别名 (Alternative Names)	NEGR1; IGLON4; Neuronal growth regulator 1; IgLON family member 4
基因 ID (Gene ID)	257194.0
蛋白 ID (SwissProt ID)	Q7Z3B1.The antiserum was produced against synthesized peptide derived from human NEGR1. AA range:172-221

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

产品名称: NEGR1 Rabbit Polyclonal Antibody
产品货号: AP Rab14544




function: May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain., similarity: Belongs to the immunoglobulin superfamily. IgLON family., similarity: Contains 3 Ig-like C2-type (immunoglobulin-like) domains., function: May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain., similarity: Belongs to the immunoglobulin superfamily. IgLON family., similarity: Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,


研究领域 (Research Area)

Cell adhesion molecules (CAMs);

图片 (Image Data)




Western blot analysis of lysates from COLO cells, using NEGR1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using NEGR1 antibody.

产品名称: NEGR1 Rabbit Polyclonal Antibody
产品货号: APRab14544



Western Blot analysis of various cells using NEGR1 Polyclonal Antibody

注意事项 (Note)

For research use only .