

产品名称: NEIL3 Rabbit Polyclonal Antibody
产品货号: APRab14547



产品概述 (Summary)

产品名称 (Production Name)	NEIL3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEIL3
别名 (Alternative Names)	NEIL3; Endonuclease 8-like 3; DNA glycosylase FPG2; DNA glycosylase/AP lyase Neil3; Endonuclease VIII-like 3; Nei-like protein 3
基因 ID (Gene ID)	55247.0
蛋白 ID (SwissProt ID)	Q8TAT5.The antiserum was produced against synthesized peptide derived from human NEIL3. AA range:491-540

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	35kDa

产品名称: NEIL3 Rabbit Polyclonal Antibody
产品货号: AP Rab14547




研究背景 (Background)

NEIL3 belongs to a class of DNA glycosylases homologous to the bacterial Fpg/Nei family. These glycosylases initiate the first step in base excision repair by cleaving bases damaged by reactive oxygen species and introducing a DNA strand break via the associated lyase reaction (Bandaru et al., 2002 [PubMed 12509226]). [supplied by OMIM, Mar 2008], function: Reports about DNA glycosylase activity are contradictory. The protein lacks a proline residue at the N-terminus that is the active site residue found in the other members of the FPG family., similarity: Belongs to the FPG family., similarity: Contains 1 FPG-type zinc finger., similarity: Contains 1 RanBP2-type zinc finger., tissue specificity: Detected in thymus and testis.,


研究领域 (Research Area)

Base excision repair;

图片 (Image Data)



Western blot analysis of lysates from HeLa, HepG2, HUVEC, and COLO cells, using NEIL3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using NEIL3 Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

产品名称: NEIL3 Rabbit Polyclonal Antibody
产品货号: APRab14547



注意事项 (Note)

For research use only .