

产品名称: Nek5 Rabbit Polyclonal Antibody
产品货号: AP Rab14554



产品概述 (Summary)

产品名称 (Production Name)	Nek5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEK5
别名 (Alternative Names)	NEK5; Serine/threonine-protein kinase Nek5; Never in mitosis A-related kinase 5; NimA-related protein kinase 5
基因 ID (Gene ID)	341676.0
蛋白 ID (SwissProt ID)	Q6P3R8.The antiserum was produced against synthesized peptide derived from human NEK5. AA range:160-209

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: Nek5 Rabbit Polyclonal Antibody
产品货号: APRab14554




研究背景 (Background)


catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Magnesium.,similarity:Belongs to the protein kinase superfamily. NEK Ser/Thr protein kinase family. NIMA subfamily.,similarity:Contains 1 protein kinase domain.,catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Magnesium.,similarity:Belongs to the protein kinase superfamily. NEK Ser/Thr protein kinase family. NIMA subfamily.,similarity:Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NEK5 Antibody. The lane on the right is blocked with the NEK5 peptide.



Western Blot analysis of various cells using Nek5 Polyclonal Antibody

注意事项 (Note)

For research use only .