

产品名称: NET-5 Rabbit Polyclonal Antibody
产品货号: APRab14574



产品概述 (Summary)

产品名称 (Production Name)	NET-5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSPAN9
别名 (Alternative Names)	TSPAN9; NET5; Tetraspanin-9; Tspan-9; Tetraspan NET-5
基因 ID (Gene ID)	10867.0
蛋白 ID (SwissProt ID)	O75954.The antiserum was produced against synthesized peptide derived from human NET-5. AA range:72-121

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	27kDa

研究背景 (Background)


产品名称: NET-5 Rabbit Polyclonal Antibody
产品货号: AP Rab14574




The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. Alternatively spliced transcripts encoding the same protein have been identified. [provided by RefSeq, Nov 2009],PTM:Glycosylated.,similarity:Belongs to the tetraspanin (TM4SF) family.,subcellular location:Colocalizes with GP6 in tetraspanin microdomains on the platelet surface.,subunit:Found in a complex with GP6.,tissue specificity:Expressed in megakaryocytes and platelets (at protein level),

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from A549 cells, using NET-5 antibody.



Western Blot analysis of various cells using NET-5 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .