

产品名称: Neurexophilin-3 Rabbit Polyclonal Antibody
产品货号: AP Rab14599



产品概述 (Summary)

产品名称 (Production Name)	Neurexophilin-3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NXPH3
别名 (Alternative Names)	NXPH3; KIAA1159; NPH3; Neurexophilin-3
基因 ID (Gene ID)	11248.0
蛋白 ID (SwissProt ID)	O95157.The antiserum was produced against synthesized peptide derived from human NXPH3. AA range:158-207

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	30kDa

研究背景 (Background)


产品名称: Neurexophilin-3 Rabbit Polyclonal Antibody
产品货号: AP Rab14599




function: May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins., PTM: May be proteolytically processed at the boundary between the N-terminal non-conserved and the central conserved domain in neuron-like cells., similarity: Belongs to the neurexophilin family., tissue specificity: Highest level in brain., function: May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins., PTM: May be proteolytically processed at the boundary between the N-terminal non-conserved and the central conserved domain in neuron-like cells., similarity: Belongs to the neurexophilin family., tissue specificity: Highest level in brain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HT-29 cells, using NXPH3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using NXPH3 antibody.

产品名称: Neurexophilin-3 Rabbit Polyclonal Antibody
产品货号: APRab14599



Western Blot analysis of various cells using Neurexophilin-3 Polyclonal Antibody

注意事项 (Note)

For research use only .