

产品名称: Neurocan Rabbit Polyclonal Antibody
产品货号: AP Rab14604



产品概述 (Summary)

产品名称 (Production Name)	Neurocan Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCAN
别名 (Alternative Names)	NCAN; CSPG3; NEUR; Neurocan core protein; Chondroitin sulfate proteoglycan 3
基因 ID (Gene ID)	1463.0
蛋白 ID (SwissProt ID)	O14594.Synthesized peptide derived from Neurocan . at AA range: 370-450

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	145kDa

研究背景 (Background)

产品名称: Neurocan Rabbit Polyclonal Antibody
产品货号: APRab14604




Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration. [supplied by OMIM, Jul 2002],function:May modulate neuronal adhesion and neurite growth during development by binding to neural cell adhesion molecules (NG-CAM and N-CAM). Chondroitin sulfate proteoglycan; binds to hyaluronic acid.,online information:Neurocan,similarity:Belongs to the aggrecan/versican proteoglycan family.,similarity:Contains 1 C-type lectin domain.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,similarity:Contains 1 Sushi (CCP/SCR) domain.,similarity:Contains 2 EGF-like domains.,similarity:Contains 2 Link domains.,tissue specificity:Brain.,

研究领域 (Research Area)

Cell Adhesion Proteins; ECM Proteins; Neuroscience; Development

图片 (Image Data)



Western Blot analysis of Jk cells using Neurocan Polyclonal Antibody

注意事项 (Note)

For research use only .