

产品名称: Neuromedin-S Rabbit Polyclonal Antibody
产品货号: AP Rab14612



产品概述 (Summary)

产品名称 (Production Name)	Neuromedin-S Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NMS
别名 (Alternative Names)	NMS; Neuromedin-S
基因 ID (Gene ID)	129521.0
蛋白 ID (SwissProt ID)	Q5H8A3.The antiserum was produced against synthesized peptide derived from human NMS. AA range:104-153

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)


产品名称: Neuromedin-S Rabbit Polyclonal Antibody
产品货号: APRab14612



This gene encodes a member of the neuromedin family of neuropeptides. The encoded preproprotein is proteolytically processed to generate a biologically active neuropeptide that plays a role in the regulation of circadian rhythm, anorexigenic action, antidiuretic action, cardiovascular function and stimulation of oxytocin and vasopressin release. [provided by RefSeq, May 2016],function:Implicated in the regulation of circadian rhythms through autocrine and/or paracrine actions.,similarity:Belongs to the NmU family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat, HepG2, HUVEC, 293, and MCF-7 cells, using NMS Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Neuromedin-S Polyclonal Antibody diluted at 1:1000

注意事项 (Note)

For research use only .