

产品名称: Neurotrimin Rabbit Polyclonal Antibody
产品货号: AP Rab14625



产品概述 (Summary)

产品名称 (Production Name)	Neurotrimin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NTM
别名 (Alternative Names)	NTM; IGLON2; NT; Neurotrimin; hNT; IgLON family member 2
基因 ID (Gene ID)	50863.0
蛋白 ID (SwissProt ID)	Q9P121.The antiserum was produced against synthesized peptide derived from human NT. AA range:235-284

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)

产品名称: Neurotrimin Rabbit Polyclonal Antibody
产品货号: AP Rab14625




This gene encodes a member of the IgLON (LAMP, OBCAM, Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009],function:Neural cell adhesion molecule.,similarity:Belongs to the immunoglobulin superfamily. IgLON family.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,


研究领域 (Research Area)

Cell Adhesion Proteins; Membrane Proteins; Neuroscience; Neurology process; Neurogenesis

图片 (Image Data)




Western blot analysis of lysates from HUVEC, HepG2, and Jurkat cells, using NT Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using NT antibody.

产品名称: Neurotrimin Rabbit Polyclonal Antibody
产品货号: APRab14625



Western Blot analysis of various cells using Neurotrimin Polyclonal Antibody

注意事项 (Note)

For research use only .