

产品名称: NF-H Rabbit Polyclonal Antibody  
产品货号: AP Rab14648



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | NF-H Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | IHC, ICC, IF, ELISA             |
| 种属反应性 (Reactivity)     | Human, Rat, Mouse               |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                            |
|------------------------|----------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | NEFH                                                                                                                       |
| 别名 (Alternative Names) | NEFH; KIAA0845; NFH; Neurofilament heavy polypeptide; NF-H; 200 kDa neurofilament protein; Neurofilament triplet H protein |
| 基因 ID (Gene ID)        | 4744.0                                                                                                                     |
| 蛋白 ID (SwissProt ID)   | P12036. The antiserum was produced against synthesized peptide derived from human NF-H. AA range: 923-972                  |

## 产品应用 (Application)

|                          |                                                           |
|--------------------------|-----------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) |                                                           |

**产品名称: NF-H Rabbit Polyclonal Antibody**  
**产品货号: AP Rab14648**




## 研究背景 (Background)

Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and functionally maintain neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene encodes the heavy neurofilament protein. This protein is commonly used as a biomarker of neuronal damage and susceptibility to amyotrophic lateral sclerosis (ALS) has been associated with mutations in this gene. [provided by RefSeq, Oct 2008],disease:Defects in NEFH are a cause of susceptibility to amyotrophic lateral sclerosis (ALS) [MIM:105400]. ALS is a neurodegenerative disorder affecting upper and lower motor neurons, and resulting in fatal paralysis. Sensory abnormalities are absent. Death usually occurs within 2 to 5 years. The etiology is likely to be multifactorial, involving both genetic and environmental factors.,function:Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber. NF-H has an important function in mature axons that is not subserved by the two smaller NF proteins.,online information:ALS genetic mutations db,polymorphism:The number of repeats is shown to vary between 29 and 30.,PTM:Phosphorylation seems to play a major role in the functioning of the larger neurofilament polypeptides (NF-M and NF-H), the levels of phosphorylation being altered developmentally and coincident with a change in the neurofilament function.,PTM:There are a number of repeats of the tripeptide K-S-P, NFH is phosphorylated on a number of the serines in this motif. It is thought that phosphorylation of NFH results in the formation of interfilament cross bridges that are important in the maintenance of axonal caliber.,similarity:Belongs to the intermediate filament family.,

## 研究领域 (Research Area)

Amyotrophic lateral sclerosis (ALS);

## 图片 (Image Data)



Immunohistochemistry analysis of NF-H antibody in paraffin-embedded human lung carcinoma tissue.

## 注意事项 (Note)

For research use only .