

产品名称: Nicalin Rabbit Polyclonal Antibody
产品货号: APRab14693



产品概述 (Summary)

产品名称 (Production Name)	Nicalin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCLN
别名 (Alternative Names)	NCLN; Nicalin; Nicastrin-like protein
基因 ID (Gene ID)	56926.0
蛋白 ID (SwissProt ID)	Q969V3.The antiserum was produced against synthesized peptide derived from human NCLN. AA range:474-523

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	63kDa

研究背景 (Background)

产品名称: Nicalin Rabbit Polyclonal Antibody
产品货号: APRab14693





function: May antagonize Nodal signaling and subsequent organization of axial structures during mesodermal patterning., similarity: Belongs to the nicastrin family., subunit: Interacts with NOMO2. Due to the strong similarity between NOMO1, NOMO2 and NOMO3, it probably also interacts with NOMO1 and NOMO3. Participates in a large protein complex, which is not related to the gamma-secretase complex., tissue specificity: Highly expressed in pancreas and skeletal muscle and, at lower levels, in heart., function: May antagonize Nodal signaling and subsequent organization of axial structures during mesodermal patterning., similarity: Belongs to the nicastrin family., subunit: Interacts with NOMO2. Due to the strong similarity between NOMO1, NOMO2 and NOMO3, it probably also interacts with NOMO1 and NOMO3. Participates in a large protein complex, which is not related to the gamma-secretase complex., tissue specificity: Highly expressed in pancreas and skeletal muscle and, at lower levels, in heart.,

研究领域 (Research Area)

Axis formation and symmetry; Developmental Biology; Embryogenesis; Morphogens;
Signal Transduction; Growth Factors/Hormones; TGF

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100 (4°, overnight). High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from

产品名称: Nicalin Rabbit Polyclonal Antibody
产品货号: APRab14693



antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .