

产品名称: NICE4 Rabbit Polyclonal Antibody
产品货号: APRab14695



产品概述 (Summary)

产品名称 (Production Name)	NICE4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBAP2L
别名 (Alternative Names)	UBAP2L; KIAA0144; NICE4; Ubiquitin-associated protein 2-like; Protein NICE-4
基因 ID (Gene ID)	9898.0
蛋白 ID (SwissProt ID)	Q14157.The antiserum was produced against synthesized peptide derived from human UBAP2L. AA range:561-610

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	115kDa

研究背景 (Background)

产品名称: NICE4 Rabbit Polyclonal Antibody
产品货号: APRab14695




PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 UBA domain.,tissue specificity:Ubiquitous.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 UBA domain.,tissue specificity:Ubiquitous.,


研究领域 (Research Area)

Cell Biology; Other Antibodies


图片 (Image Data)



Western blot analysis of lysates from HT-29 and 293 cells, using UBAP2L Antibody. The lane on the right is blocked with the synthesized peptide.




Western blot analysis of the lysates from Jurkat cells using UBAP2L antibody.



产品名称: NICE4 Rabbit Polyclonal Antibody
产品货号: APRab14695



Western Blot analysis of NIH-3T3 K562 NIH-3T3 293T cells using NICE4 Polyclonal Antibody diluted at 1: 500



Western Blot analysis of HEPG2 using NICE4 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .