

产品名称: NKp44 Rabbit Polyclonal Antibody
产品货号: APRab14724



产品概述 (Summary)

产品名称 (Production Name)	NKp44 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCR2 NCR2; LY95; Natural cytotoxicity triggering receptor 2; Lymphocyte antigen 95
别名 (Alternative Names)	homolog; NK cell-activating receptor; Natural killer cell p44-related protein; NK-p44; NKp44; CD336
基因 ID (Gene ID)	9436.0
蛋白 ID (SwissProt ID)	O95944.The antiserum was produced against synthesized peptide derived from the Internal region of human NCR2. AA range:61-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	44kDa

产品名称: NKp44 Rabbit Polyclonal Antibody
产品货号: APRab14724




研究背景 (Background)

function:Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.,similarity:Belongs to the natural cytotoxicity receptor (NCR) family.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,subunit:Interacts with TYROBP/DAP12.,tissue specificity:Selectively expressed by activated NK cells and by in vitro cultured (i.e. activated) TCRg/d lymphoid cells.,function:Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.,similarity:Belongs to the natural cytotoxicity receptor (NCR) family.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,subunit:Interacts with TYROBP/DAP12.,tissue specificity:Selectively expressed by activated NK cells and by in vitro cultured (i.e. activated) TCRg/d lymphoid cells.,


研究领域 (Research Area)

Natural killer cell mediated cytotoxicity;

图片 (Image Data)




Western Blot analysis of PC-12 cells using NKp44 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000




Immunohistochemical analysis of paraffin-embedded human-Lymph-nodes, antibody was diluted at 1:100

产品名称: NKp44 Rabbit Polyclonal Antibody
产品货号: APRab14724



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100

注意事项 (Note)

For research use only .