

产品名称: Nkx-2.4 Rabbit Polyclonal Antibody
产品货号: APRab14729



产品概述 (Summary)

产品名称 (Production Name)	Nkx-2.4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NKX2-4
别名 (Alternative Names)	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D
基因 ID (Gene ID)	644524.0
蛋白 ID (SwissProt ID)	Q9H2Z4.The antiserum was produced against synthesized peptide derived from human NKX24. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: Nkx-2.4 Rabbit Polyclonal Antibody
产品货号: AP Rab14729




研究背景 (Background)

function: Probable transcription factor., similarity: Belongs to the NK-2 homeobox family., similarity: Contains 1 homeobox DNA-binding domain., function: Probable transcription factor., similarity: Belongs to the NK-2 homeobox family., similarity: Contains 1 homeobox DNA-binding domain.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Transcription Factors

图片 (Image Data)



Western blot analysis of the lysates from K562 cells using NKX24 antibody.



Western Blot analysis of various cells using Nkx-2.4 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .