

产品名称: Nkx-2.6 Rabbit Polyclonal Antibody
产品货号: APRab14731



产品概述 (Summary)

产品名称 (Production Name)	Nkx-2.6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NKX2-6
别名 (Alternative Names)	NKX2-6; NKX2F; Homeobox protein Nkx-2.6; Homeobox protein NK-2 homolog F
基因 ID (Gene ID)	18092.0
蛋白 ID (SwissProt ID)	A6NCS4.The antiserum was produced against synthesized peptide derived from human NKX26. AA range:91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	28kDa

产品名称: Nkx-2.6 Rabbit Polyclonal Antibody
产品货号: APRab14731





研究背景 (Background)


This gene encodes a homeobox-containing protein that belongs to the NK-2 homeobox family. This protein is a vertebrate homolog of *Drosophila* homeobox-containing protein called 'tinman', which has been shown to be essential for development of the heart-like dorsal vessel. In conjunction with related gene, NKX2-5, this gene may play a role in both pharyngeal and cardiac embryonic development. Mutations in this gene are associated with persistent truncus arteriosus. [provided by RefSeq, Aug 2011],

研究领域 (Research Area)

图片 (Image Data)



产品名称: Nkx-2.6 Rabbit Polyclonal Antibody
产品货号: APRab14731



Western Blot analysis of various cells using Nkx-2.6 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .